

PDE7A Positive Control for STJ502369 peptide (STJ504667)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

PRODUCT PROPERTIES

Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	18583
Gene Symbol	Pde7a
UniProt ID	PDE7A_MOUSE
Specificity	This is positive control is recommended for use in combination with PDE7A antibody STJ502369.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---