

PDE9A Positive Control for STJ502390 peptide (STJ504679)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	WB
PRODUCT PROPERTIES	
Dilution Range	WB: Apply 15-20 μ l per lane
Formulation	Provided as 100 μ L ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	5152
Gene Symbol	PDE9A
UniProt ID	PDE9A_HUMAN
Specificity	This is positive control is recommended for use in combination with PDE9A antibody STJ502390.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.