

## FGFR4 Positive Control peptide (Phospho) (STJ504692)

### GENERAL INFORMATION

|              |              |
|--------------|--------------|
| Product Type | Biomolecules |
| Applications | WB           |

### PRODUCT PROPERTIES

|                     |                                                                 |
|---------------------|-----------------------------------------------------------------|
| Dilution Range      | WB: 1:500                                                       |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">2264</a>        |
| Gene Symbol | <a href="#">FGFR4</a>       |
| UniProt ID  | <a href="#">FGFR4_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|