

## Smad1 Positive Control peptide (Phospho) (STJ504748)

### GENERAL INFORMATION

Product Type Biomolecules

Applications WB

### PRODUCT PROPERTIES

Dilution Range WB: 1:500-1:1, 000

Storage Instruction Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.