

## MAP3K12 Positive Control peptide (Phospho) (STJ504750)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

### PRODUCT PROPERTIES

Dilution Range	WB: 1:500
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">7786</a>
Gene Symbol	<a href="#">MAP3K12</a>
UniProt ID	<a href="#">M3K12_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---