

PDE3A Blocking Peptide peptide (Phospho) (STJ504775)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
---------------------	---

TARGET INFORMATION

Gene ID	54611
Gene Symbol	Pde3a
UniProt ID	PDE3A_MOUSE
Immunogen Region	200-300

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---