

PDE5A Positive Control peptide (Phospho) (STJ504787)

GENERAL INFORMATION

Product Type **Biomolecules**

Applications **WB**

PRODUCT PROPERTIES

Dilution Range **WB: 1:500**

Storage Instruction **Store at -20°C for long term storage. Avoid freeze-thaw cycles.**

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**