

## Src Blocking Peptide peptide (Phospho) (STJ504822)

### GENERAL INFORMATION

|              |                                   |
|--------------|-----------------------------------|
| Product Type | Biomolecules                      |
| Applications | Immunodepletion/Immunocompetition |

### PRODUCT PROPERTIES

|                     |                                                                 |
|---------------------|-----------------------------------------------------------------|
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
|---------------------|-----------------------------------------------------------------|

### TARGET INFORMATION

|                  |                           |
|------------------|---------------------------|
| Gene ID          | <a href="#">6714</a>      |
| Gene Symbol      | <a href="#">SRC</a>       |
| UniProt ID       | <a href="#">SRC_HUMAN</a> |
| Immunogen Region | 400-500                   |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|