

ROR1 Blocking Peptide for STJ502825 peptide (STJ504869)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Immunogen Region	900-918
Specificity	This blocking peptide is recommended for use in combination with ROR1 antibody, STJ502825
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.