

Adenylate Cyclase 10 Blocking Peptide for STJ502868 peptide (STJ504883)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	271639
Gene Symbol	Adcy10
UniProt ID	ADCYA_MOUSE
Immunogen Region	200-250
Specificity	This blocking peptide is recommended for use in combination with Adenylate Cyclase 10 antibody, STJ502868
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.