

Sall4C Blocking Peptide for STJ502872 peptide (STJ504884)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	ELISA/IHC/IP/WB

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-50
Specificity	This blocking peptide is recommended for use in combination with Sall4C antibody, STJ502872

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---