

TIA1 Blocking Peptide for STJ503256 peptide (STJ504993)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Dilution Range	WB: IHC: 1:50 ICC: 1:50 IF: 1:50 CFM: 1:50
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7072
Gene Symbol	TIA1
UniProt ID	TIA1_HUMAN
Immunogen Region	414-464
Specificity	This blocking peptide is recommended for use in combination with TIA1 antibody, STJ503256

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---