

TLR1 Blocking Peptide for STJ503276 peptide (STJ505003)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	7096
Gene Symbol	TLR1
UniProt ID	TLR1_HUMAN
Immunogen Region	370-420
Specificity	This blocking peptide is recommended for use in combination with TLR1 antibody, STJ503276
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.