

TM6SF2 Blocking Peptide for STJ503314 peptide (STJ505013)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	DB/ELISA/IHC/IP/WB

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	107770
Gene Symbol	Tm6sf2
UniProt ID	TM6S2_MOUSE
Immunogen Region	220-270
Specificity	This blocking peptide is recommended for use in combination with TM6SF2 antibody, STJ503314

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---