

ZMPSTE24 Blocking Peptide peptide (STJ505082)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Concentration	0.55 µg/µl
Dilution Range	WB: DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10269
Gene Symbol	ZMPSTE24
UniProt ID	FACE1_HUMAN
Immunogen Region	215-245

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---