

Relaxin 1 Positive Control for STJ502816 peptide (STJ505101)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	WB
PRODUCT PROPERTIES	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	295144
Gene Symbol	Rxtp1
UniProt ID	RXFP1_RAT
Specificity	This is positive control is recommended for use in combination with Relaxin 1 antibody STJ502816.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.