

## SPHK1 Positive Control for STJ503116 peptide (STJ505170)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

### PRODUCT PROPERTIES

Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">20698</a>
Gene Symbol	<a href="#">Sphk1</a>
UniProt ID	<a href="#">SPHK1_MOUSE</a>
Specificity	This is positive control is recommended for use in combination with SPHK1 antibody STJ503116.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---