

SUR1 Positive Control for STJ503162 peptide (STJ505180)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	WB
PRODUCT PROPERTIES	
Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	25559
Gene Symbol	Abcc8
UniProt ID	ABCC8_RAT
Specificity	This is positive control is recommended for use in combination with SUR1 antibody STJ503162.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.