

VEGFR1 Blocking Peptide peptide (STJ505319)

GENERAL INFORMATION

| | |
|--------------|-----------------------------------|
| Product Type | Biomolecules |
| Applications | Immunodepletion/Immunocompetition |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 2321 |
| Gene Symbol | FLT1 |
| UniProt ID | VGFR1_HUMAN |
| Immunogen Region | 1250-1300 |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|