

MeCP2 isoform 1 Blocking Peptide for STJ501738 peptide (STJ505349)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	4204
Gene Symbol	MECP2
UniProt ID	MECP2_HUMAN
Immunogen Region	1-50
Specificity	This blocking peptide is recommended for use in combination with MeCP2 isoform 1 antibody, STJ501738
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.