

FGFR3 Positive Control peptide (Phospho) (STJ505364)

GENERAL INFORMATION

Product Type Biomolecules

Applications WB

PRODUCT PROPERTIES

Dilution Range WB: 1:500

Storage Instruction Store at -20°C for long term storage. Avoid freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.