

## hCG beta Peptide Blocking Peptide peptide (STJ505378)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

### PRODUCT PROPERTIES

Dilution Range	WB: IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:250
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">114335</a>
Gene Symbol	<a href="#">CGB1</a>
UniProt ID	<a href="#">CGB1_HUMAN</a>
Immunogen Region	142-165

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---