

STAT1 Blocking Peptide peptide (Phospho) (STJ505409)

GENERAL INFORMATION

| | |
|--------------|-----------------------------------|
| Product Type | Biomolecules |
| Applications | Immunodepletion/Immunocompetition |

PRODUCT PROPERTIES

| | |
|---------------------|-----------------------------------------------------------------|
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |
|---------------------|-----------------------------------------------------------------|

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 6772 |
| Gene Symbol | STAT1 |
| UniProt ID | STAT1_HUMAN |
| Immunogen Region | 700-750 |

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|