

ZFP36L1 Blocking Peptide for STJ503573 peptide (STJ505474)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	677
Gene Symbol	ZFP36L1
UniProt ID	TISB_HUMAN
Immunogen Region	200-250
Specificity	This blocking peptide is recommended for use in combination with ZFP36L1 antibody, STJ503573

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---