

## IL18RAP Positive Control for STJ501437 peptide (STJ505483)

### GENERAL INFORMATION

|              |              |
|--------------|--------------|
| Product Type | Biomolecules |
| Applications | WB           |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Dilution Range      | WB: 1:500   |
| Formulation         | Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.   |

### TARGET INFORMATION

|             |   |
|-------------|---|
| Gene ID     | <a href="#">8807</a>  |
| Gene Symbol | <a href="#">IL18RAP</a>   |
| UniProt ID  | <a href="#">I18RA_HUMAN</a>   |
| Specificity | This is positive control is recommended for use in combination with IL18RAP antibody STJ501437. |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|