

BMPR1A Blocking Peptide for STJ500281 peptide (STJ505490)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Dilution Range	WB: IHC: 1:50-1:150
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	657
Gene Symbol	BMPR1A
UniProt ID	BMR1A_HUMAN
Immunogen Region	460-510
Specificity	This blocking peptide is recommended for use in combination with BMPR1A antibody, STJ500281

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---