

Syncytin 2 Positive Control for STJ503177 peptide (STJ505534)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

PRODUCT PROPERTIES

Dilution Range	WB: 1:500
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	405754
Gene Symbol	ERVFRD-1
UniProt ID	SYCY2_HUMAN
Specificity	This is positive control is recommended for use in combination with Syncytin 2 antibody STJ503177.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---