

Cell Surface Glycoprotein MUC18 Blocking Peptide for STJ501820 peptide (STJ505575)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	78967
Gene Symbol	Mcam
UniProt ID	MUC18_RAT
Immunogen Region	330-353
Specificity	This blocking peptide is recommended for use in combination with Cell Surface Glycoprotein MUC18 antibody, STJ501820

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---