

HIPK1 Blocking Peptide for STJ501341 peptide (STJ505581)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Concentration	0.5-1.5 $\mu\text{g}/\mu\text{l}$
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	204851
Gene Symbol	HIPK1
UniProt ID	HIPK1_HUMAN
Immunogen Region	700-750
Specificity	This blocking peptide is recommended for use in combination with HIPK1 antibody, STJ501341
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.