

## GRK6 Positive Control for STJ501288 peptide (STJ505602)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

### PRODUCT PROPERTIES

Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemelli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2870</a>
Gene Symbol	<a href="#">GRK6</a>
UniProt ID	<a href="#">GRK6_HUMAN</a>
Specificity	This is positive control is recommended for use in combination with GRK6 antibody STJ501288.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---