

## AMPK alpha 1 Blocking Peptide for STJ500075 peptide (STJ505610)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

### PRODUCT PROPERTIES

Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">65248</a>
Gene Symbol	<a href="#">Prkaa1</a>
UniProt ID	<a href="#">AAPK1_RAT</a>
Immunogen Region	500-550
Specificity	This blocking peptide is recommended for use in combination with AMPK alpha 1 antibody, STJ500075

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---