

GPS2 Blocking Peptide for STJ501285 peptide (STJ505643)

GENERAL INFORMATION	
Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition
PRODUCT PROPERTIES	
Dilution Range	WB: 1:2, 500-1:5, 000
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	2874
Gene Symbol	GPS2
UniProt ID	GPS2_HUMAN
Immunogen Region	250-300
Specificity	This blocking peptide is recommended for use in combination with GPS2 antibody, STJ501285
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.