

## DNER Blocking Peptide for STJ500779 peptide (STJ505647)

| GENERAL INFORMATION    |                                                                                                                                                                |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Biomolecules                                                                                                                                                   |
| Applications           | Immunodepletion/Immunocompetition                                                                                                                              |
| PRODUCT PROPERTIES     |                                                                                                                                                                |
| Dilution Range         | WB: 1:2, 500-1:5, 000                                                                                                                                          |
| Formulation            | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction    | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                |
| TARGET INFORMATION     |                                                                                                                                                                |
| Gene ID                | <a href="#">92737</a>                                                                                                                                          |
| Gene Symbol            | <a href="#">DNER</a>                                                                                                                                           |
| UniProt ID             | <a href="#">DNER_HUMAN</a>                                                                                                                                     |
| Immunogen Region       | 700-737                                                                                                                                                        |
| Specificity            | This blocking peptide is recommended for use in combination with DNER antibody, STJ500779                                                                      |
| ADDITIONAL INFORMATION |                                                                                                                                                                |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                        |