

OR14I1 Positive Control for STJ505969 peptide (STJ505745)**GENERAL INFORMATION**

Product Type	Biomolecules
Applications	WB

PRODUCT PROPERTIES

Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	401994
Gene Symbol	OR14I1
UniProt ID	Q14I1_HUMAN
Specificity	This is positive control is recommended for use in combination with OR14I1 antibody STJ505969.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---