

CRALBP Blocking Peptide for STJ505978 peptide (STJ505750)

GENERAL INFORMATION

| | |
|--------------|-----------------------------------|
| Product Type | Biomolecules |
| Applications | Immunodepletion/Immunocompetition |

PRODUCT PROPERTIES

| | |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Formulation | Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice. |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|---------------------------------------------------------------------------------------------|
| Gene ID | 6017 |
| Gene Symbol | RLBP1 |
| UniProt ID | RLBP1_HUMAN |
| Immunogen Region | 260-317 |
| Specificity | This blocking peptide is recommended for use in combination with CRALBP antibody, STJ505978 |

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|