

## Androgen Receptor Blocking Peptide peptide (Phospho) (STJ505895)

### GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

### PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
---------------------	---

### TARGET INFORMATION

Gene ID	<a href="#">367</a>
Gene Symbol	<a href="#">AR</a>
UniProt ID	<a href="#">ANDR_HUMAN</a>
Immunogen Region	550-600

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---