

Anti-PDE2A Antibody antibody (N-Term) (STJ505918)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CM/ELISA/ICC/IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.7 µg/µl
Conjugation	Unconjugated
Dilution Range	WB: 1:250-1:1, 000 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:250-1:1, 000 ICC: 1:250-1:1, 000 IF: 1:150-1:1, 000 CFM: 1:250-1:1, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, HEPES, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	N-Term
------------------	--------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---