

## Anti-Rabbit Anti-Chicken IgG HRP antibody (STJ505919)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CM/ELISA/IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Chicken

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---