

## Anti-TUFM antibody (325-375) (STJ506047)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/IP/WB        |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Polyclonal  |
| Concentration       | 0.5-1.5 $\mu\text{g}/\mu\text{l}$   |
| Conjugation         | Unconjugated  |
| Purification        | Affinity Purified   |
| Dilution Range      | WB: 1:500<br>ELISA: 1:50, 000<br>IP: 1:200  |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Heps, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG   |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.   |

### TARGET INFORMATION

|                  |                            |
|------------------|----------------------------|
| Gene ID          | <a href="#">7284</a>       |
| Gene Symbol      | <a href="#">TUFM</a>       |
| UniProt ID       | <a href="#">EFTU_HUMAN</a> |
| Immunogen Region | 325-375                    |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|