

Anti-CGRP antibody (350-410) {Biotin} (STJ506113)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/IP/WB |
| Host / Source | Rabbit |
| Reactivity | Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality | Polyclonal |
| Concentration | 0.5-1.5 $\mu\text{g}/\mu\text{l}$ |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500-1:5,000 ELISA: 1:50,000 IP: 1:200 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|-----------------------------|
| Gene ID | 12310 |
| Gene Symbol | Calca |
| UniProt ID | CALCA_MOUSE |
| Immunogen Region | 350-410 |

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|