

Anti-TESK2 antibody (200-250) (STJ506191)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/IP/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.5-1.5 Mu g/Mu l
	Conjugation	Unconjugated
	Purification	Affinity Purified
	Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:200
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	10420
	Gene Symbol	TESK2
	UniProt ID	TESK2_HUMAN
	Immunogen Region	200-250
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.