

TNF alpha Positive Control for STJ506269 peptide (STJ506271)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

PRODUCT PROPERTIES

Formulation	Provided as 100 uL ready-to-use, in SDS-PAGE sample buffer (Laemmli's buffer) containing Tris, pH 6.8, 1 % SDS, Glycerol and Bromophenolblue blue as tracking dye. The sample is reduced by adding 2% beta mercaptoethanol. The protein concentration is
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7124
Gene Symbol	TNF
UniProt ID	TNFA_HUMAN
Immunogen Region	100-150
Specificity	This is positive control is recommended for use in combination with TNF alpha antibody STJ506269.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---