

Anti-TMEM262 antibody (1-75 aa) {Biotin} (STJ506272)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/IP/WB
	Host / Source	Rabbit
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.5-1.5 Mu g/Mu l
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	ELISA 1:10,000 ICC 1:50-1:100 IF 1:50-1:100 IHC 1:50-1:100 IP 1:200 WB 1:500
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	433215
	Gene Symbol	Tmem262
	UniProt ID	TM262_MOUSE
	Immunogen Region	1-75 aa
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.