

Anti-SARS CoV-2 Spike Glycoprotein Brazil variant-Brazil variant antibody (400-450) (STJ506298)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Rabbit
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	ELISA 1:10, 000 ICC 1:50-1:100 IF 1:50-1:100 IHC 1:50-1:100 IP 1:200 WB 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, HEPES, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	400-450
------------------	---------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---