

Anti-Phospho-PLK4 antibody (300-350) (STJ506388)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/IP/WB |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Polyclonal |
| Concentration | 0.5-1.5 Mu g/Mu l |
| Conjugation | Unconjugated |
| Purification | Affinity Purified |
| Dilution Range | ELISA 1:10, 000 IP 1:200 WB 1:500 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, HEPES, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|----------------------------|
| Gene ID | 10733 |
| Gene Symbol | PLK4 |
| UniProt ID | PLK4_HUMAN |
| Immunogen Region | 300-350 |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|