

ENAH Blocking Peptide for STJ506413 peptide (Phospho) (STJ506416)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	Immunodepletion/Immunocompetition

PRODUCT PROPERTIES

Dilution Range	ELISA 1:10,000 IP 1:200 WB 1:500
Formulation	Liquid form at 2.5mg/ml concentration in PBS. Up to 5% DMSO can be added. Orders with >1mg can be supplied in lyophilized powder form, or in buffer of choice.
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13800
Gene Symbol	Enah
UniProt ID	ENAH_MOUSE
Immunogen Region	710-760
Specificity	This blocking peptide is recommended for use in combination with ENAH antibody, STJ506413

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---