

ENAH Positive Control peptide (Phospho) (STJ506417)

GENERAL INFORMATION

Product Type	Biomolecules
Applications	WB

PRODUCT PROPERTIES

Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
---------------------	---

TARGET INFORMATION

Gene ID	13800
Gene Symbol	Enah
UniProt ID	ENAH_MOUSE
Immunogen Region	710-760

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---