

## Anti-ENAH/MENA antibody (710-760) (STJ506418)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Rabbit
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	ELISA 1:10, 000 IP 1:200 WB 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">13800</a>
Gene Symbol	<a href="#">Enah</a>
UniProt ID	<a href="#">ENAH_MOUSE</a>
Immunogen Region	710-760

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---