

Anti-ADCY2 antibody (C-Term) {Biotin} (STJ506424)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	0.5-1.5 Mu g/Mu l
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	ELISA 1:10, 000 IP 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene Symbol	Adcy2
UniProt ID	ADCY2_MOUSE
Immunogen Region	C-Term
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.