

Anti-SKOR1 antibody (N-Term) {Biotin} (STJ507767)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/IP/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Conjugation	Biotin
	Purification	Affinity Purified
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	24473
	Gene Symbol	Htr1a
	UniProt ID	5HT1A_RAT
	Immunogen Region	N-Term
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.